



Liquid property control system for LVAD ISO 5198 testing and mock circulatory loop simulations.

Jack R. Friedberg

PI: Charles E. Taylor, Ph.D.



CAHL

Cajun Artificial Heart Laboratory




UNIVERSITY
OF
LOUISIANA
L a f a y e t t e

Outline


- What was the Liquid Property Control Unit (LPCU) designed for?
- Design and structure of the unit
- Methods and materials for the operation of the unit
- Overview of the program that governs the operations of the LPCU
- Methods and materials used for collecting the results
- Results
- Conclusions
- Future work

Background

- Mock circulatory loops (MCLs) test viability of LVADs
- Blood viscosity and density varies so LVADs must be tested at a range of viscosities and densities to ensure safety
- Previously, the liquid in MCLs were changed by hand
- LPCU allows for round-the-clock scheduling of simulations and simplistic means of documenting changes



LPCU Design Overview



Refutas Equation

$$VBN = 14.534 \ln[\ln(v + 0.8)] + 10.975$$

$$VBN_{Blend} = (X_A \times VBN_A) + \dots + (X_N \times VBN_N)$$



$$v = \exp \left[\exp \left(\frac{VBN_{Blend} - 10.975}{14.534} \right) \right] - 0.8$$

$v = \text{Viscosity}$



$VBN = \text{Viscosity Blending Number}$

Materials and Methods for Construction of LPCU

- A Graphical User Interface (GUI) made with MathWorks Matlab sends user input to the Arduino





- An Arduino Mega 2560 with a DFRobot motor shield computes the liquid to be dispensed and controls the pumps



Materials and Methods

- Fluids are pumped into and out of the loop with Bartendro peristaltic pumps
- A relay system activates solenoid valves that allow for all three base liquids to be drawn through the same pump
- The three base fluids stored in plastic pales:
 - 60:40 glycerin-water mixture by weight
 - Pure distilled water
 - Saturated NaCl mixture




Design of the LPCU


Pumps and valves



Arduino with motor shield



Relay system



Base liquids



Graphical User Interface

Working_Gui

Fluid Property Control

| Properties | | Status | |
|---------------------|----------------------|--|----------------------|
| Existing Viscosity: | <input type="text"/> | <input type="text"/> cSt | <input type="text"/> |
| Desired Viscosity: | <input type="text"/> | <input type="text"/> cSt | <input type="text"/> |
| Desired Density: | <input type="text"/> | <input type="text"/> Kg/m ³ | <input type="text"/> |
| Total Volume: | <input type="text"/> | <input type="text"/> L | <input type="text"/> |

COM Port #

Connect

No Connection

Disconnect

Run

Working_Gui

Fluid Property Control

| Properties | | Status | |
|---------------------|---------------------------|--|----------------------|
| Existing Viscosity: | <input type="text"/> 4 | <input type="text"/> cSt | <input type="text"/> |
| Desired Viscosity: | <input type="text"/> 5 | <input type="text"/> cSt | <input type="text"/> |
| Desired Density: | <input type="text"/> 1125 | <input type="text"/> Kg/m ³ | <input type="text"/> |
| Total Volume: | <input type="text"/> 8 | <input type="text"/> L | <input type="text"/> |

COM Port #

8





Connect

Connected


Disconnect

Run

How the LPCU is Used




Arduino Code Flow Diagram



Materials and Methods for Verifying Results

- Samples of the resultant fluid are taken using a volumetric pipette and stored in glass vials
- A balance is used to verify the density of the sample
- A pocket refractometer is used to verify the refractive index of the sample
- A viscometer is used to verify the viscosity of the sample



Results

- Using the instruments on the previous page, confirmed the properties of 3 known mixtures
- Density and viscosity of glycerin/water solutions
 - 40/60 glycerin/water
 - 50/50
 - 60/40

| Glycerin/Water Solutions | Theoretical Density | Measured Density | Theoretical Viscosity | Measured Viscosity |
|--------------------------|---------------------|------------------|-----------------------|--------------------|
| 40/60 | 1.11 | 1.12 g/mL | 3.12 cSt | 3.14 cSt |
| 50/50 | 1.14 | 1.16 g/mL | 4.98 cSt | 4.94 cSt |
| 60/40 | 1.26 g/mL | 1.25 g/mL | 8.72 cSt | 7.9 cSt* |

Results

- 3 containers with a glycerin water mix were made
- It was determined that the LPCU accurately dispenses the volume of desired fluids as measured in a graduated container

| | Initial | | Target | | Final | |
|--------------------|------------|-----------|------------|-----------|------------|-----------|
| | Density | Viscosity | Density | Viscosity | Density | Viscosity |
| Container 1 | 1.102 g/mL | 2.81 cSt | 1.125 g/mL | 3.50 cSt | 1.114 g/mL | 3.52 cSt |
| Container 2 | 1.104 g/mL | 3.00 cSt | 1.125 g/mL | 4.00 cSt | 1.252 g/mL | 3.86 cSt |
| Container 3 | 1.103 g/mL | 3.17 cSt | 1.125 g/mL | 4.00 cSt | 1.119 g/mL | 3.83 cSt |

Discussion

- The Arduino is correctly calculating the fluid to be drained and added
- Arduino accurately dispenses the amount of each fluid that the Arduino requests
- The accuracy of the base fluids is the main concern at this point
 - Slow separation of water/glycerin mixture
 - Unreliable density of NaCl mixture from temperature changes
 - Mold developing in liquids over time

Conclusion

- The LPCU is well on it's way to being fully operational and only minor adjustments need to be made to the base fluids
 - Constant stirring
 - Working to keep temperature constant
 - Using UV lights to kill bacteria
- Designed to be easily customizable and can be used for other applications



Future Direction

Short Term Goals

- Further work with base fluids to ensure accuracy of the final mixture
- Improving the GUI to better inform the user on the status of the Unit and document liquid changes

Long Term Goals

- Shear thinning



CAHL
Cajun Artificial Heart Laboratory

Thank You

Jack R. Friedberg

Research Assistant, Mechanical Engineering

jackfriedberg@gmail.com